

Applications to Medical Imaging from diagnostic to therapy



By P. Le Du

dapnia

cea

saclay

patrick.le-du@cea.fr




Protontherapy (PSI)

X-Rays, the fastest technology transfer example



- On November 8, 1895 Röntgen discovered X-Rays
- On November 22, 1895 he takes the first image of his wife's hand




Röntgen received the first Nobel prize in physics in 1901

Medical Imaging: a multidisciplinary approach

Physics



Mathematics



Medicine




Chemistry



Biology



Informatics




Medical Imaging Modalities

X-ray radiology

X-ray Radiography

Computed Tomography (CT)


Tomo-Densitometry (DXA)



Emission tomography

Single Photon Emission Computerized Tomography (SPECT)

Positron Emission Tomography (PET)




Magnetic Resonance (MR)

Magnetic Resonance Imaging (MRI)

Magnetic Resonance Spectroscopy (MRS)


Functionnal MRI (fMRI)



Ultrasonic imaging



The various types of imaging



New!

- Complementary !
- Depends on what you want to see

NMR & PET Images of Epilepsy

NMR




PET



- NMR "Sees" Structure with 0.5 mm Resolution
- PET "Sees" Metabolism with 5.0 mm Resolution but with very high sensitivity (picomolar level)

Positron Emission Tomography principle

- Functional imaging
- Molecular tracers with doped beta + emitters
 - ^{11}C (20 min), ^{15}O (2 min), ^{13}N (10 min), ^{18}F (2 h)
 - Produced by a 18 MeV Proton cyclotron
 - The most common $\rightarrow ^{18}\text{F} \Rightarrow ^{18}\text{FDG}$ fluoro-deoxy-glucose
 - Sign the degree of activity of an organ hungry of glucose
- annihilation positron with an electron
 - emission of two 511 keV photons back to back



Physical constraints & limits

Emission



β^+ flight path: ^{18}F : $E = 650 \text{ keV} \Rightarrow \langle \lambda \rangle = 2 \text{ mm}$

accolinearity $\gamma\text{--}\gamma$: kinematic : $\langle \Theta \rangle \text{few mRad}$


Detection

Detector resolution


Bias reconstruction



Source of errors in the source




Some history




Jeavons, Townsend & al.

			PET III 1975
			ECAT II 1977
			NeuroECAT 1978
			ECAT 931 1985
			ECAT EXACT HR ⁺ 1995



Today



Clinical PET Imaging


- Whole body PET PET $\rightarrow \Phi=90\text{cm}$ - FOV 20 cm
 - Oncology \rightarrow market increase 30% per year !




NaI curved
CPET (Philips)




CTI
Siemens
LSO



Detection block
Crystals $4 \times 4 \times 20$ (or 30) mm^3
Block 8×8 crystals , 2×2 PM's



Patient treated for a colon cancer
and revealing under PET/CT scan
an additional breast cancer



C-PET
Philips
1997



PET -CT Biograph

2006




GEMS Discovery

Medical PET examples : Lung Cancer

New Cases 185,800/yr
Deaths 163,700/yr




HOW DOES PET COMPARE?



(Gould et al., JAMA, 2001)

Others types of common cancers

Colon



New Cases 147,500/yr
Deaths 57,100/yr


Breast



New Cases 220,900/yr
Deaths 28,900/yr

New Cases 212,600/yr
Deaths 40,200/yr


Prostate



New Cases 220,900/yr
Deaths 28,900/yr


PET for brain diseases

Alzheimer



The PET shows decreased metabolism early in the disease!

Parkinson



The PET scan showed abnormal glucose metabolism in the back of the right hemisphere. Following surgical removal of the dysfunctional brain area, the child was seizure free.

■ MicroPET → $\Phi=20\text{cm}$ - FOV few cm


- Radio pharmacology
- Tracer development




LabPET



ClearPET




APD PET



RatCap



31 g mouse
1 mCi ^{18}F -





Whole-body
FDG-PET scan
250 g rat
(Sherbrooke APD)




Dedicated PET



Mammography (CLEARPEM)



Prostate PET



On Line PET for hadrontherapy



μ PET vs whole body PET

→ different requirements

■ High Spatial resolution

→ fundamental

- Objective ~ 1mm or less
- Today → 1,2 mm

■ High sensitivity

- Less Compton event
- Small dose

■ Parallax correction

→ Depth Of Interaction Technique



■ High Efficiency (>85%)

■ Good Spatial Resolution (<5 mm)

■ Low Cost (<\$100/cm²)

■ Short Dead Time (<1 μ s)


■ High Timing Resolution (<5 ns fwhm)

■ Good Energy Resolution (<100 keV fwhm)




Why PET ?

Similarities and differences




Calorimeter

HEP



$M_{\text{Higgs}} = 100 \text{ GeV}$



PET Camera


Biomedical Imaging

Similarities

Geometry and granularity
 Detector (Crystals & scintillator)
 Sensor (PM, APD)
 Electronics: Fast (40 MHz), compact
 Event rate & Data volume (Gbit/s)

Differences

Energy range (10 GeV-511 keV)
 No synchronisation
 --> free running electronics
 Multiple vertices








From HEP to Medical

Where **techniques** are transferred to developments in bio- medical field
Medical Imaging has only partially benefited from new technologies developed for telecommunications and High Energy Physics detectors

- **New scintillating crystals and detection materials** →
 - CMS (WPbO₄) → Luap ...(Crystal Clear col)
- **Photodetectors : Highly segmented and compact** → PMT → APD → SiPM
 - APD : SSC/SDC (1991) → CMS (1996) → MicroTEP → TEP
- **Electronics & signal treatment** → Highly integrated
 - Fast, low noise, low power preamp
 - Digital filtering and signal analysis
- **Trigger/DAQ** →
 - High level of parallelism and event filtering algorithms
 - Pipeline and parallel read-out, trigger and on-line treatment
- **Computing**
 - Modern and modular simulation software using worldwide recognized standards (GEANT)

From HEP detectors to Medical Imaging

Requirements for HEP EM calorimetry


1. **Crystals** Technology transfer 
 - High density ($> 6 \text{ g/cm}^3$)
 - Fast emission ($< 100 \text{ ns}$), visible spectrum
 - Moderate to high light yield
 - High radiation resistance
2. **Photodetectors** Technology transfer 
 - Compact
 - High quantum efficiency and high gain
 - High stability
3. **Readout electronics** Technology transfer 
 - Fast shaping, low noise
 - Highly integrated
4. **Intelligent and parallel DAQ** Technology transfer 
 - Reduce dead time
5. **Software** Technology transfer 
 - Accurate Monte Carlo simulation
6. **General design** Technology transfer 
 - Compact integration of a large number of channels ($> 10'000$)

Requirements for Medical Imaging

1. **Crystals**
 - High density ($> 7 \text{ g/cm}^3$)
 - Fast emission ($< 100 \text{ ns}$), visible spectrum
 - High light yield
 - Moderate radiation resistance
2. **Photodetectors**
 - Compact
 - High quantum efficiency and high gain
 - High stability
3. **Readout electronics**
 - Fast shaping, low noise
 - Highly integrated
4. **Intelligent and parallel DAQ**
 - Reduce dead time
5. **Software**
 - Accurate Monte Carlo simulation
6. **General design**
 - Compact integration of a large number of channels ($> 10'000$)

The Present: Combine anatomic and functional informations

morphology **metabolism**



David Townsend
CERN: 1970-78
Université de Genève

UPSM Pittsburgh
and

Ronald Nutt
(CTS – CTI)



Today state of the art: TEP + CT

Spiral CT


Survey

CT

University of Pittsburgh
PET/CT scanner




scatter correction
attenuation correction




TEP

FUSION





Merged image
TEP/CT₂₂

In 60% of case the
bimodalimages have a
consequence in the
future therapeutic
treatment protocol!...






Patient Radiation Dose is Limited!




- Image Noise Is Limited by Counting Statistics
- Cannot Increase too much Source Strength

Volumetric CT



< 0,4 sec/ rotation
Organ in a sec (17 cm/sec)
Whole body < 10 sec


20 to 50 mSv
Standard radiography
0.1 mSv




R & D examples

1- Crystals

CMS PbWO₄ production



Crystal Clear LuAP production




LuYAP
Crystals

Scintillators for PET


	1962	1977	1995	1999	2001	LaBr ₃
	NaI	BGO	GSO:Ce	LSO:Ce	LuAP:Ce	
Density (g/cm ³)	3.67	7.13	6.71	7.40	8.34	
Atomic number	51	75	59	66	65	
Photofraction	0.17	0.35	0.25	0.32	0.30	
Decay time (ns)	230	300	30-60	35-45	17	
Light output (hv/MeV)	43000	8200	12500	27000	11400	
Peak emission (nm)	415	480	430	420	365	
Refraction index	1.85	2.15	1.85	1.82	1.97	
Hygroscopic	Yes	No	No	No	No	
Natural radioactivity	No	No	Yes	Yes	Yes	

No Scintillator with Superior Properties in *All Aspects*


Detectors → crystals




NaI curved
CPET (Philips)



Detection block
Crystals 4x4x20 (or 30) mm³
Block 8 x 8 crystals , 2 x 2 PM's



GEMS
BGO (Bicron)




ADAC Philips
GSO



CTI Siemens
LSO

2- New pixellised Photodetectors

CMS



Hamamatsu single channel APD


LHCB



HPD tube
manufactured at
CERN: 2048 channels

BrainPET

ClearPEM



Hamamatsu 32
channels APD array

Opera



Hamamatsu H7546
64 channel PMT

ClearPET




Mammography




Hamamatsu
PM flat panel

State of the art technologies


A Small Animal PET device
Licensed to industry (ClearPET)




(courtesy of P.Lecoq/CERN)



Hamamatsu 32 channels
APD array





The Rat Conscious Animal PET scanner, J.F. Pratte, et al., BNL




Ring containing 12 block detectors
Up to two layers of 2x2 x5 mm deep
LSO crystals with APDs and integrated
readout electronics

APD
(Hamamatsu
S8550)

LSO array




Counting rate estimate




From
D. Crossetto

4- Pipeline Architectures

LHC



Future PET



Digitisation


Pipeline

Event builder

Dead Time

Time during which the detector can not accept and record new events

Dead time sources:



$$R = \text{Rate (s}^{-1}\text{)}$$


$$\delta = \text{Absolute dead time per event (s)}$$

$$DT = \text{Relative dead time (\%)}$$

$$AR = \text{Acquisition rate (s}^{-1}\text{)}$$

$$DT = 1 - e^{-R\delta}$$


$$AR = R \cdot e^{-R\delta}$$




Pipeline Architectures

- Architectures based on multi-event storage and pipeline trigger processing
- Allow to minimize the data acquisition dead time

Example with 4 channels and trigger $A+B+C+D > \text{threshold}$



Proposal of Front End Architecture




- ◆ Trigger logic processes "raw fast information"
- ◆ Free-running sampling ADC
- ◆ Digital filter used to extract pulse amplitude and high resolution timing
- ◆ Pipelined processing architecture to avoid deadtimes
- ◆ Only one "channel" to compute either the energy and time



Towards an innovative read out electronics concept (INNOTEP)

- ◆ Truly pixelising the detector yields to a considerable number of channels to be considered, each having its own reading electronics
- ◆ size and speed have been a real issue for years but progress in the microelectronics field have made ASICs of high integration readily available
- ◆ besides this, they appear as a cost effective solution
- ◆ no possible CFD implementation on chip (shared constant network / derivation are noisy...!) in agreement with the expected time resolution.
- ◆ High resolution TDCs require complex architecture, large surface area and appear no
→ need to find another solution for time measurement ... to be inspired by HEP experiments

Next → Time-of-Flight in PET ?




- Can localize source along line of flight.
- Time of flight information reduces **noise** in images.
- Time of flight cameras built in the 80's with BaF2 & CsF.
- These scintillators forced compromises that prevented TOF from flourishing
- Today new crystals (LaBr3) and new MHz electronics/DAQ






■ **Objective : < 500 ps Timing Resolution**

Whole-Body TOF Simulations

2×10^6 Trues, 1×10^6 Randoms, Attenuation Included
OP-OSEM w/ TOF Extensions, 2 Iterations, 14 Subsets



(1:2:3 body:liver:tumor)




Clear Improvement Visually

*Data courtesy of Mike Casey, CPS Innovations

The future of medical imaging


- Faster exams
- Movement correction
 - Breathing
 - Cardiac beating
 - Digestive bolus
- Dynamics
- Quantification
- True multimodality
- Reduce dose to patient



IMPROVE

- Spatial resolution
- Timing resolution
- Sensitivity
- Signal/Noise ratio

Molecular Imaging in Medicine & Biology



■ At present

- Tissue sample during surgery or after biopsy
- Molecular analysis to define the tumor characteristics
→ molecular signature

■ In the near future

- Tissue sample
by imaging before surgery
→ molecular signature


PET Molecular Imaging




Accelerator

◆ ^{15}O
◆ ^{13}N
◆ ^{11}C
◆ ^{18}F


2-[F-18]Fluoro-2-Deoxy-D-Glucose (FDG)



Healthcare
Solutions



Pre-Clinical



Clinical

Future Challenges

Whole-Body High-Resolution 'Minute' PET

- Very large number of channels (20 → 80 cm FOV)
~ 300 k channels (2x2 mm² pixels)
- High trigger rate ~ 10 MHz
(10 mCi, 20% sensitivity)
- High data rate
~ 10 Gbyte/s (1 kbyte event size)
- Large number of events
~ 10⁹ events (10⁶ voxels, 1000 events/voxel)
- Large data volume per image
~ 1000 Gbytes (list mode)
- High computer power for image reconstruction

Improve spatial resolution ?

- Near theoretical limit = few mm (4mm with FDG)
- Can Increase SNR by Reducing Backgrounds
- Keep exam time short (30 min → few min!)



Sensitivity → significant room for improvement!

- Compact and hermetic design → large Field Of View
- Fast and high light yield crystal (LaBr₃)
- Fast and low noise electronics with TOF capability
- Built-in intelligence in the Data Acquisition system
- Make the best use of "good" events (TRIGGER)
- Use Compton events instead of rejecting them?
- Efficiently throw away "bad" events (better timing resolution!)

5- Simulation

Higgs event at LHC (CMS) with Geant4



ClearPET with GATE: Geant4 Application
for Tomographic Emission



GATE : Geant4 Application for Tomographic Emission


Monte-Carlo simulation allowing to :

- ✓ define geometries
(size, materials,...)
- ✓ define sources
(geometry, nature, activity)
- ✓ choice of physical process
(low energy package of G4)
- ✓ follow track point by point




GATE specificities:

- ✓ CERN GEANT4 libraries
- ✓ Time modelling
(sources, movement, random...)
- ✓ Script language (avoid C++)
- ✓ Code interactivity
- ✓ Sharing development






Simulation of decaying sources




^{15}O (2 min)

^{11}C (20 min)



Radiotherapy :Context



• Radiotherapy:

- Use ionizing radiation to treat and cure tumors
- Frequent treatment (2/3 of cases) → 160 000 /year in France.
- Efficient treatment → 30 to 40% of recovery




→ Failures due to radioresistant tumors!

- Allow good quality of life → tolerance
 - non invasive, itinerant and without important physical effects.
- Cheap → 5% of the cancer budget (500 M€)
- Essentially X (Linear accelerators) & photons (curietherapy)
- Why Radiotherapy X is NOT 100 % efficient?
 - Local irradiation → 100 Gy = 90 % of sterilization
 - Complication < 5 %
 - Tolerance of saine tissue is the limiting factor of Rx
 - Technological progress = improve ratio: **tumor dose / sain tissue dose**

Hadrontherapy accelerators the rationale


200 MeV - 1 nA
protons

4800 MeV – 0.1 nA
carbon ions
which can control
radioresistant
tumours




From U. Amaldi

How to irradiate the tumor ?




Spread Bragg peak



■ Treatment in depth → combine

- Energy modulation → Scan the energy to make a Spread Out Bragg Peak (SOBP) that spans the tumor
- Intensity modulation

Comparison with Photons



state-of-the-art photon therapy
(IMRT)



Proton therapy
(spot scanning)

Carbon Ion Radiotherapy

73M Lt. Nasal Cavity Malignant Melanoma T4N0M0 57.6GyE/16fr/4w



Before



2 months after RT

Specific Aims of PET in Heavy Ion Therapy

- **Verification of treatment plan *in vivo* prior to dose delivery**
 - ≈ measure the beam stopping point with high spatial accuracy (<1 mm)
 - *High spatial resolution*
 - *High sensitivity*
- **Monitor dose distribution “in-beam”**
 - ≈ measure b^+ activity with high quantitative accuracy
 - *High sensitivity (→ high SNR !)*
 - *High count rate capability*
 - *Fast image reconstruction (real time..?)*

R & D → Some areas of interest

■ Beam production

- New synchrotrons!
- Active scanning of the tumor

■ Beam monitoring and dosimetry →

- Solid state and micropattern detectors, fast electronics ...


■ Patient dosimetry

- In beam PET → On line monitoring of deposited doses need high resolution and high sensitive PET in real time environment


■ Proton CT →

- Particle telescope (solid state or gaseous) and calorimeters, fast electronics

GSI treatment planning using « active » scanning





Treatment room
at GSI




In Situ PET Control of light Ion therapy physical basis


- Activation of the irradiated tissue by positron emitters (^{11}C , ^{10}C , ^{15}O ...)
 - By product of the therapeutic irradiations
- Detection by PET
 - In beam, in situ, non invasive
- Relation between dose and β^+ activities
 - Beam position, particle range
- Cross sections
 - $^{12}\text{C} + ^{16}\text{O} \rightarrow ^{15}\text{O}$: 84 mb
 - $^{12}\text{C} + ^{12}\text{C} \rightarrow ^{11}\text{C}$: 56 mb
 - $^{12}\text{C} + ^{12}\text{C} \rightarrow ^{13}\text{C}$: 5 mb




Clinical results at GSI



Treatment plan



Predicted β^+ -activity





Measured β^+ -activity


- Range verification for portal with high penetration depth and inhomogeneous target volume

Clinical results at GSI


- Range verification in delicate situations:
 - Stopping of the beam in front of organ at risk



Treatment plan



Predicted
 β^+ -activity



Measured
 β^+ -activity

Thanks to

- R. Lecomte (CHU Sherbrooke)
- P. Lecoq (CERN)
- R. Ferrand (CP Orsay)
- Prof. J.N. Talbot (Hôpital Tenon- Paris)
- J. Remilieux & M. Bajard - (ETOILE project)
- Prof J.P. Gérard (Nice)
- U. Amaldi (Pavia project)
- W. Enghardt & K. Parodi (FZR-GSI)
- M. Grossmann (PSI)
- W.W. Moses (LBL)
- C.Woody (BNL)
- R.Aymar (CERN)and many others